# EndoPilot<sup>2</sup>

## Processing instructions (cleaning, disinfection and sterilisation)

Attachment for the user manual

<b>Contents</b>		
1.	Processing instructions for thermostable products .....	4
1.1.	General Information	4
1.2.	Pre-treatment	4
1.3.	Automated cleaning / disinfection (with RDG (thermal washer disinfectors))	5
1.4.	Manual cleaning and disinfection	5
1.5.	Inspection	6
1.6.	Maintenance / Using service oil	6
1.7.	Packaging	6
1.8.	Sterilisation	6
1.9.	Steam sterilisation	6
1.10.	Storage	6
1.11.	Reusability	6
2.	Processing instructions for thermolabile products.....	8
2.1.	General principles	8
2.2.	Pre-treatment	8
2.3.	Manual cleaning and disinfection	8
2.4.	Inspection	9
2.5.	Maintenance	9
2.6.	Packaging and storage	9
2.7.	Sterilisation	9
2.8.	Reusability	9

Illustration 1

Illustration 2

Illustration 3

Illustration 4



Illustration 4e

Contact wire, housing, spring, threaded bush, knob, contact

You can dismantle the patented file clamp (4e) for reprocessing.  
A universal wrench for ultrasonic-instruments can be used as a tool for assembly (with a wrench size of 3.2 mm).

Illustration 5

Illustration 6

Illustration 7

### Description of the single parts and assignment of the relevant processes

ill.	Ref.	Description	Chapter
1	110 2010 110 2011	Only the control unit with a touch screen	Ch. 2
2	110 2203	Power supply with primary plug Input: 100 – 240 V AC Output: 12 V DC 1.50 A	Ch. 2
3	109 2361	Bluetooth wireless foot switch, single pedal	Ch. 2
4	109 2311 v06	Apex cable set	Ch. 2
	109 2312 109 2314 109 2318 109 2315 109 2316	4a – Measuring cable with plug 4b – Lip clip 4c – Cap for the plug socket 4d – Cable for file clamp 4e – File clamp (can be dismantled)	Ch. 1 Ch. 1 Ch. 1 Ch. 1
	110 2303	The file clamp can be dismantled (see illustration 4e) To dismantle, the contact is unscrewed and removed from the knob. You can clean all parts individually. At the end of assembly, firmly tighten the contact again. Caution. Check if the device is functioning properly. Loose parts may fall out and enter the patient's mouth.  4f – Retainer for apex cable (mounted on the device)	Ch. 1
5	109 0126	Contra-angle handpiece with apex measurement / fully insulated 1:1 gear, with ISO-E coupling	A3*
6	109 0112	Motor with apex measuring contact, LED power indicator and ISO-E connection	Ch. 2
7	109 0151	DownPack handpiece with a LED display	Ch. 2
	540 5173 364 2901	7a – Screw cap 7b – blue O-ring	Ch. 1 Ch. 1

A3\* Note the user manual and processing instructions for the contra-angle:

See document 609 4110

# 1. Processing instructions for thermostable products

## 1.1. General Information

Read the processing instructions carefully and completely. Follow the instructions section by section.

All products/product components must be cleaned, disinfected and sterilised before each application; this applies in particular to the first use after delivery, since all products/product components are delivered non-sterile (cleaning and disinfection take place after removal of the protective transport packaging; sterilisation after repackaging). Effective cleaning and disinfection are indispensable prerequisites for effective sterilisation.

As part of your responsibility for the sterility of the products/product components during application, please note that

- only sufficiently device-specific and product-specific validated methods are used for cleaning, disinfection and sterilisation,
- the equipment used (RDG `cleaning and disinfection device´ and sterilizer) is regularly maintained and correctly checked and
- the validated parameters are adhered to during each cycle.

Material resistance:

When selecting cleaning agents and disinfectants, please make sure that the following components are not included:

- organic, mineral and oxidising acids (minimum permissible pH value 5.5)
- strong alkalis (maximum permissible pH value 12, slightly alkaline cleaner recommended)
- organic solvents (e.g. alcohols, ethers, ketones, benzines)
- Oxidising agents (e.g. hydrogen peroxides)
- Halogens (chlorine, iodine, bromine)
- aromatic/halogenated hydrocarbons

Always refrain from cleaning products/product components with metal brushes or steel wool.

All products/product components may only be exposed to temperatures not exceeding 138°C (280°F)!

If possible, an automated method (RDG (cleaning and disinfection device)) should be used for cleaning and disinfection. Due to its significantly lower effectiveness and reproducibility, a manual method, even when using an ultrasonic bath, should only be used if an automated method is not available.

Pre-treatment must be carried out in both cases.

Please also observe the legal regulations valid in your country and the hygiene regulations of your medical practice or hospital. This applies in particular to the different specifications regarding effective prion inactivation.

Additional or deviating specifications apply to some products/product components (see chapter `Special instructions´ in the attachment).

## 1.2. Pre-treatment

Immediately after the application (within a maximum of 2 h) all visible residuals must be removed from the products/product components:

Care must be taken when selecting the combined cleaning- and disinfection-agent to ensure

- that these are suitable for the cleaning and disinfection of instruments made of metals and plastics in principle,
- that these are suitable for ultrasonic cleaning (no foam development),
- that these have a tested effectiveness (e.g. VAH/DGHM listing or CE marking) and
- that these are compatible with the products/product components (see chapter `Material resistance´, section 1.1.).

The concentrations, temperatures, residence times and rinsing specifications specified by the manufacturer of the cleaning/ disinfection- agent must be adhered to in any case. Use only freshly prepared solutions, only DI water (The bacterial load must correspond at least to that of drinking water with a maximum of 100 CFU. The endotoxin load should be low). Only use a soft, clean and lint-free cloth and/or filtered air for drying.

**Procedure:**

1. Remove externally visible dirt with a disposable cloth dipped in the cleaning and disinfection bath.
2. Dismantle the products/product components as far as possible (see also the chapter `Special instructions´).
3. Place the products/product components in the cleaning/disinfection solution in the ultrasonic bath (temperature <35°C/95°F) and activate the ultrasound for the specified residence time (but not less than 5 min).
4. After switching off the ultrasound, manually remove all visible dirt in a separate cleaning bath. Use a clean soft brush (1 minute brushing) that you only use for this purpose; never use metal brushes or steel wool. (For recommendations on the brush type, see `Special instructions´ in the attachment)
5. Rinse or flush thoroughly with DI water for at least 1 minute (The bacterial load must correspond at least to that of drinking water with a maximum of 100 CFU and with a low endotoxin load at the same time).

**Now continue to section 1.3. If an automated method is not available continue to 1.4**

### 1.3. Automated cleaning / disinfection (with RDG (thermal washer disinfectant))

Care must be taken when selecting the (RDG) thermo washer disinfectant to ensure

- that the thermo washer disinfectant has a tested effectiveness (CE marking) in principle,
- that a tested program for thermal disinfection ( $A_0$  value > 3000 or – for older devices – at least 5 min at 90°C/194°F) is used wherever possible (risk of disinfectant residues on the products/product components in the event of chemical disinfection),
- that the program used is suitable for the products/product components and contains sufficient rinsing cycles,
- that only DI water is used for rinsing (The bacterial load must correspond at least to that of drinking water with a maximum of 100 CFU and with a low endotoxin load at the same time) (The DI supply of the device must be checked regularly, as bacteria can accumulate in cartridge systems. Filter systems to reduce the bacterial load are recommended),
- that the air used for drying is filtered (oil-free, low bacterial count and low particle count) and
- that the thermo washer disinfectant is regularly maintained and correctly checked.

Care must be taken when selecting the cleaning agent system used to ensure

- that these are suitable for the cleaning of instruments made of metals and plastics in principle,
- that – if no thermal disinfection is used – a suitable disinfectant with tested effectiveness (e.g. VAH/DGHM-listing or CE marking) is additionally used and that this is compatible with the cleaning agent used and
- that the chemicals used are compatible with the products/product components (see chapter ‘Material resistance’ in section 1.1.).

The concentrations, temperatures, residence times and rinsing specifications specified by the manufacturer of the cleaning agent and disinfectant, must be adhered to in any case.

#### Procedure:

1. Place the dismantled products/product components (see chapter ‘Pre-treatment’) in the thermo washer disinfectant. Make sure that the products/product components do not touch each other in the process. Use a “small parts basket” for smaller components (see ‘Special instructions’).
2. Start the program. e.g. Vario-TD (Miele)
3. Remove the products/product components from the thermo washer disinfectant at the end of the program.
4. Inspect and pack the products/product components as soon as possible after removal (see the chapters ‘Inspection’, ‘Maintenance’ and ‘Packaging’, if necessary after additional drying in a clean place).

*An independent accredited and recognised (§ 15 (5) MPG) test laboratory (Medical Device Services – DR: ROSSBERGER GmbH, Gilching) has provided proof of the general suitability of the products/product components for effective automated cleaning and disinfection by using the RDG G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the cleaning agent Neodisher MediClean forte (Dr. Weigert GmbH & Co. KG, Hamburg). The selected parameters corresponded to the Vario-TD program from Miele. The method described above was considered and the combined cleaning/disinfection agent ID 213 (DÜRR Dental, Bietigheim-Bissingen) with a concentration of 2% was used for preliminary cleaning.*

### 1.4. Manual cleaning and disinfection

Care must be taken when selecting the combined cleaning- and disinfection-agent to ensure

- that these are suitable for the cleaning and disinfection of instruments made of metals and plastics in principle,
- that these are suitable for ultrasonic cleaning (no foam development),
- that these have a tested effectiveness (e.g. VAH/DGHM listing or CE marking) and
- that these are compatible with the products/product components (see chapter ‘Material resistance’ in section 1.1).

The concentrations, temperatures, residence times and rinsing specifications specified by the manufacturer of the cleaning/ disinfection- agent must be adhered to in any case. Use only freshly prepared solutions, only DI water (The bacterial load must correspond at least to that of drinking water with a maximum of 100 CFU. The endotoxin load should be low). Only use a soft, clean and lint-free cloth and/or filtered air for drying.

#### Procedure

##### Cleaning:

1. Place the dismantled products/product components (see chapter ‘Pre-treatment’) in the ultrasonic bath filled with the cleaning/ disinfection-agent for the specified residence time so that the products/product components are sufficiently covered. Make sure that the products/product components do not touch each other in the process and activate the ultrasound (do not use any heating function) for the specified residence time (but not less than 5 min).
2. After switching off the ultrasound, support cleaning in a separate cleaning bath by completely brushing off all internal and external surfaces with a soft brush. (For recommendations on the brush type, see ‘Special instructions’ in the attachment)
3. Then remove the products/product components from the cleaning bath and rinse or flush these thoroughly with DI water at least three times (for at least 1 minute each) (The bacterial load must correspond at least to that of drinking water with a maximum of 100 CFU and with a low endotoxin load at the same time).
4. Inspect the products/product components (see the chapters ‘Inspection’ and ‘Maintenance’).

##### Disinfection

5. Place the dismantled, cleaned and inspected products/product components in the ultrasonic bath in a fresh cleaning/disinfectant solution for the specified residence time so that the products/product components are sufficiently covered. Make sure that the products/product components do not touch each other in the process.
6. Then remove the products/product components from the disinfection bath and rinse or flush these thoroughly with DI water at least five times (for at least 1 minute each) (The bacterial load must correspond at least to that of drinking water with a maximum of 100 CFU and with a low endotoxin load at the same time).
7. Dry the products/product components by blowing off/blowing out with sterile filtered compressed air.
8. Pack the products/product components as soon as possible after drying (see the chapter ‘Packaging’, if necessary after additional drying in a clean place).

*An independent accredited and recognised (§ 15 (5) MPG) test laboratory (Medical Device Services – DR: ROSSBERGER GmbH, Gilching) has provided proof of the general suitability of the products/product components for effective manual cleaning and disinfection by using the combined cleaning/disinfection agent ID 213 (DÜRR Dental, Bietigheim-Bissingen) for both the cleaning and disinfection steps. The method described above was considered and the combined cleaning/disinfection agent ID 213 (DÜRR Dental, Bietigheim-Bissingen) with a concentration of 2% was used for preliminary cleaning. The Sonorex Super RK 514 H from Bandelin was used as the ultrasonic device.*

## 1.5. Inspection

After cleaning or cleaning/disinfection, check all products/product components for corrosion, damaged surfaces, chipping, dirt and discolouration, and remove damaged products/product components (for the numerical restriction on reuse see the chapter `Reusability`). Products/product components that are still dirty must be cleaned and disinfected again.

## 1.6. Maintenance / Using service oil

Reassemble the dismantled product if required. (see the chapter `Special instructions`).

Instrument oils must not be used with this product. (see the chapter `Special instructions`).

## 1.7. Packaging

Please pack the products/product components in disposable sterilisation packaging (one-way packaging) which corresponds to the following requirements (material/process):

- DIN EN ISO 11607
- Suitable for steam sterilisation (temperature resistance up to min. 138°C (280°F) sufficient vapour permeability)
- Products/product components or sterilisation packaging are sufficiently protected against mechanical damages

## 1.8. Sterilisation

Only the sterilization methods listed below are to be used for sterilising the products/product components; other sterilisation methods are not permitted.

## 1.9. Steam sterilisation

- Fractionated vacuum method<sup>1</sup> or gravitation method<sup>2</sup> (with sufficient product drying<sup>3</sup>)
- Steam steriliser according to DIN EN 13060/DIN EN 285
- Validated according to DIN EN ISO 17665 (valid IQ/OQ (picking) and product-specific performance appraisal (PQ))
- Maximum sterilisation temperature 134°C (273°F; plus tolerance according to DIN EN ISO 17665)
- Sterilisation time (exposure time at sterilisation temperature):

Country	Fractionated vacuum method	Gravitation method
Germany/Europe	Minimum of 5 min <sup>4</sup> at 132°C (270°F)/134°C (273°F), minimum drying time of 20 min <sup>3</sup>	not recommended

<sup>1</sup> At least three vacuum steps

<sup>2</sup> Use of the less effective gravitation method is only permitted if the fractionated vacuum method is not available. It requires significantly longer sterilisation times and must be validated by the user for specific products, devices, methods and parameters.

<sup>3</sup> The actual drying time required directly depends on parameters that are the sole responsibility of the user (load configuration and density, steriliser condition,...) and must therefore be determined by the user. Nevertheless, drying times should not be less than 20 minutes.

<sup>4</sup> or 18 min (prion inactivation, not relevant for USA)

Flash sterilisation is generally not permitted.

Do not use hot air sterilisation, radiation sterilisation, formaldehyde or ethylene oxide sterilisation, or plasma sterilisation.

*An independent accredited and recognised (§ 15 (5) MPG) test laboratory (Medical Device Services – DR: ROSSBERGER GmbH, Gilching) has provided proof of the general suitability of the products/product components for effective steam sterilisation by using the steam steriliser HST 6x6x6 (Zirbus technology GmbH, Bad Grund) and using the fractionated vacuum method and the special oil spray item no. 00-108-005 (Dr. Schumacher GmbH, Melsungen. Typical conditions in clinics and doctor's practices as well as the method described above were considered here.*

## 1.10. Storage

After sterilisation, the products/product components must be stored in the sterilisation packaging in a dry and dust-free manner.

## 1.11. Reusability

The products/product components can – with due care and provided that they are undamaged and free of dirt – be reused at least up to the number indicated in the chapter `Special instructions`. The user is responsible for any further use beyond this. The use of damaged and/or dirty products/product components is not permitted. All liability is excluded in the event of non-observance.

Special instructions (here you will find further information on individual parts that require a special procedure)

Item no.	Item name	Brush	Special/additional procedure for				Packaging	Sterilisation	The maximum permissible number of cycles	The corresponding recommendation classification KRINKO/RK/BfA/M recommendation (Germany only, if used as intended)
			Pre-treatment	Manual cleaning/ disinfection	Automated cleaning/ disinfection	Maintenance/ Assembly				
109 2316	File clamp EndoPilot	Standard and conical interdental brush of type IDB-GK, green, 3-6.5 mm, item no. 266668 ORBIS Dental, Münster  for thin diameters in the housing and the threaded bush: cylindrical interdental brush of type IDB-W, white, 2.8 mm, item no. 266666 ORBIS Dental, Münster	Remove the cable, unscrew the knob and separate the parts, also push the threaded bush out of the knob (if necessary, use a dental probe)  then: immerse, brush inside and outside (small holes with interdental brushes, thin diameter in housing and threaded bush with the fine, cylindrical type IDB-W 2.8 mm),  Ultrasound	completely dismantled: immerse, brush outside and inside (small inside holes with interdental brushes, thin diameter in housing and threaded bush with the fine, cylindrical type IDB-W 2.8 mm),  Ultrasound	completely dismantled: "small parts basket"	Assembly: 1. Insert the threaded bush into the knob, loosely screw in the contact 2. Insert the contact wire into the housing and press the housing onto the table in a vertical position 3. Place the spring and knob onto the contact wire standing in a vertical position 4. Align the knob, press it down and tighten the contact firmly. (wrench with a width over flats of 3.2 mm) Also see illustration 4e and table: 'Description of the single parts' row 4	Pre-assembled, transparent sterilisation packaging (simple)	Pre-assembled, standard procedure	300	Critical A/B
109 2315	File cable (for file clamp)	Standard, and conical interdental brush of type: IDB-GK, green, 3-6.5 mm, item no. 266668 ORBIS Dental, Münster	Immerse into the cleaning solution, brush, (inner surfaces of the socket with conical interdental brush)  Ultrasound	brush (inside holes with conical interdental brush),  ultrasound	in the "small parts basket", openings facing downwards	Do not oil	Transparent sterilisation packaging (simple)	Standard procedure	300	Critical A/B
109 0151	Down-Pack handpiece parts:  Screw cap 540 5173 and O-ring 364 2901	Standard, and conical interdental brush of type: IDB-GK, green, 3-6.5 mm see above	Dismantled, immerse into the cleaning solution, brush, (inner surfaces of the cap with conical interdental brush)  Ultrasound	immerse dismantled, brush inside and outside (small inside holes with conical interdental brush),  ultrasound	Dismantled, "small parts basket",	Do not oil  Change the O-rings after 50 cycles	not pre-assembled  Transparent sterilisation packaging (simple)	not pre-assembled  Standard procedure	300 (O-rings cycles)  50	Critical A/B

## 2. Processing instructions for thermolabile products

### 2.1. General principles

All products/product components must be cleaned and disinfected before each application; this also applies in particular to the first use after delivery, as all products/product components are delivered non-sterile (cleaning and disinfection after removal of the protective transport packaging). Effective cleaning and disinfection are indispensable prerequisites for safe use.

As part of your responsibility for the safety of the products/product components during application, please note that only sufficient product-specific validated methods for cleaning and disinfection are used and that the validated parameters are adhered to for each cycle.

Material resistance:

When selecting disinfectants, please make sure that the following components are not included:

- Acids or alkalis (disinfectants for residue-free surface disinfection on an alcoholic basis are recommended)
- other organic solvents (e.g. aldehydes, acetone, ethers, benzines)
- Oxidising agents (e.g. hydrogen peroxides)
- Halogens (chlorine, iodine, bromine)
- aromatic/halogenated hydrocarbons
- Oils (including instrument oils)

Always refrain from cleaning products/product components with metal brushes or steel wool.

All products/product components may only be exposed to temperatures not exceeding 60°C (140°F)!

Please also observe the legal regulations valid in your country and the hygiene regulations of your medical practice or hospital. Additional or deviating specifications apply to some products/product components (see the chapter 'Special instructions').

### 2.2. Pre-treatment

Immediately after the application (within a maximum of 2 h) all visible residuals must be removed from the products/product components and the products/product components must be disinfected:

**Procedure:**

**Pre-treatment**

1. Separate components that require separate cleaning, disinfection and sterilisation according to the specific processing instructions (see 'Special instructions' in the attachment).

### 2.3. Manual cleaning and disinfection

Selecting the disinfectant

Care must be taken when selecting the disinfectant to ensure

- that these are suitable for the residue-free disinfection of medical surfaces made of metals and plastics in principle (alcoholic basis, without fragrances and perfumes),
- that these have a tested effectiveness (e.g. VAH/DGHM listing or CE marking)
- that the chemicals used are compatible with the products/product components (see the chapter 'Material resistance' in section 2.1.).

The concentrations, temperatures and residence times specified by the manufacturer of the disinfectant must be adhered in any case. Only use a fresh solution and a soft, clean and lint-free disposable cloth for drying.

**Procedure:**

**Cleaning**

1. Remove the remaining products/product components from all accessible places in the relevant areas (see the chapter 'Special instructions' in the attachment) with a disposable wipe soaked in at least 50 ml of disinfectant.
2. Repeat step 1 until the surface appears free of visible residues.
3. Inspect the products/product components (see the chapters 'Inspection' and 'Maintenance').

**Disinfection**

4. Wet the products/product components on all accessible places in the relevant areas (see the chapter 'Special instructions' in the attachment) with a disposable wipe soaked in at least 50 ml of disinfectant and allow the disinfectant to have an effect according to the residence time specified by the manufacturer of the disinfectant (but not less than 5 minutes, moisten again if it dries too quickly).
5. Repeat step 4.
6. Allow the surface to completely dry.

*An independent officially accredited and recognised (§ 15 (5) MPG) test laboratory (Medical Device Services – DR: ROSSBERGER GmbH, Gilching) has provided proof of the general suitability of the products/product components for effective pre-treatment/manual cleaning and disinfection by using the disinfectant Bacillol AF (Paul Hartmann AG, Heidenheim). The method described above was considered here.*

## 2.4. Inspection

After cleaning or cleaning/disinfection, check all products/product components for cracks, damaged surfaces, dirt and discolouration and remove damaged products/product components (for the numerical restriction on reuse see the chapter `Reusability`). Products/product components that are still dirty must be cleaned and disinfected again.

## 2.5. Maintenance

Instrument oils must not be used.

## 2.6. Packaging and storage

Packaging is not required.

Repeat the cleaning/disinfection (see the chapter `Cleaning/Disinfection`) at the beginning of each day of treatment.

## 2.7. Sterilisation

Sterilisation is generally not permitted.

## 2.8. Reusability

The products/product components can – with due care and provided they are undamaged and free of dirt – be reused at least up to the number indicated in the chapter `Special instructions`. The user is responsible for any further use or the use of damaged and/or dirty products/product components. All liability is excluded in the event of non-observance.

**Special instructions** (here you will find further information on individual parts that require a special procedure)

Item no.	Item name	Relevance/frequency (according to the risk-based approach)	Special/additional procedure for		Sterilisation	The maximum permissible number of cycles	The corresponding recommendation for classification KRINKO/RKI/BfA/M (Germany only, if used as intended)
			Pre-treatment, manual cleaning and disinfection	Maintenance/ Assembly			
109 2312	EndoPilot apex connection cable (from the EndoPilot apex cable set, item no. 109 2311)	Mandatory after each use: - All surfaces Not permitted: - Plug/socket (interior)	Standard	Do not oil	Sterilisation not permitted	1000	Non-critical
109 2315	EndoPilot file cable (from the EndoPilot apex cable set, item no. 109 2311)	<b>This component is thermostable and must be steam sterilised</b> (see specific processing instructions for steam sterilisation, see section 1) Separate cleaning, disinfection and sterilisation (no wipe disinfection is intended)	Standard	Do not oil	Sterilisation not permitted	1000	No patient contact
109 2316	EndoPilot apex file clamp (from the EndoPilot apex cable set, item no. 109 2311)	<b>This component is thermostable and must be steam sterilised</b> (see specific processing instructions for steam sterilisation, see section 1) Disassemble the file clip and then perform separate cleaning, disinfection and sterilisation (no wipe disinfection is intended)					
109 2314	EndoPilot apex lip clip (from the EndoPilot apex cable set, item no. 109 2311)	<b>This component is thermostable and must be steam sterilised</b> (see specific processing instructions for steam sterilisation, see section 1) Separate cleaning, disinfection and sterilisation (no wipe disinfection is intended)					
109 2318	EndoPilot cap for the plug socket (from the EndoPilot apex cable set, item no. 109 2311)	<b>This component is thermostable and must be steam sterilised</b> (see specific processing instructions for steam sterilisation, see section 1) Separate cleaning, disinfection and sterilisation (no wipe disinfection is intended)					
110 2010 or 110 2011	EndoPilot <sup>2</sup> control box (and extension modules, when mounted)	Mandatory after each use: Touch screen For additional safety (recommended at least once per working day): - All other surfaces Not permitted: - Plug/socket (interior)	Standard	Do not oil	Sterilisation permitted	1000	No patient contact
110 2301 110 2302	Holder for handpieces right side Holder for handpieces left side		Standard Not permitted				

Special instructions

Item no.	Item name	Relevance/frequency (according to the risk-based approach)	Special/additional procedure for		Sterilisation	The maximum permissible number of cycles	The corresponding recommendation for KRINKO/RKI/BfArM (Germany only, if used as intended)
			Pre-treatment, manual cleaning and disinfection	Maintenance/ Assembly			
109 0112	Endo-Pilot Endo motor	Mandatory after each use: - before dismantling the contra-angle handpiece to be processed separately: Contra-angle handpiece and front part (until the separating line) of the motor - after dismantling the contra-angle handpiece to be processed separately: Front part (until the separating line) of the motor  For additional safety (recommended at least once per working day): Front (after dismantling the contra-angle handpiece to be processed separately) - All other surfaces Not permitted: - Plug/socket (interior)	Standard	Do not oil, Connect the contra-angle directly before the usage	Sterilisation not permitted	1000	Non-critical/semi-critical A
			Standard				
			Not permitted				
109 0151	D-Pack handpiece	Mandatory after each use: - before dismantling the heating needle to be processed separately and the front nut: Heating needle/front nut and front part (until the separating line) of the handpiece - after dismantling the heating needle to be processed separately and the front nut with O-ring: Front part (until the separating line) of the handpiece  For additional safety (recommended at least once per working day): Front (after dismantling the heating needle to be processed separately) - All other surfaces Not permitted: - Plug/socket (interior)	Standard	Do not oil, connect the heating tip directly before the usage	Sterilisation not permitted	1000	Non-critical/semi-critical A
			Standard				
			Not permitted				
		<b>These components are thermostable and must be steam sterilised (see specific processing instructions for steam sterilisation, see section 1)</b> Separate cleaning, disinfection and sterilisation (no wipe disinfection is intended)					

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Schlumbohm GmbH & Co. KG